



LINC

**IN.PACT DCB results from
Real world experience of a
Single vascular surgeon
Prospective registry**

Kishore Sieunarine
Vascular Surgeon
Perth WA



LINC

SFA/Pop Treatment Options

TASC A to D

POBA Initial success, restenosis up to 80% 12m

BMS: Absolute, Astron, Fast, Fact, Resilient, 4Ever 75% at 12m

DURABILITY 200, VIASTAR 65% 12M

SUPERB (SUPERA) Calcification 86% 12M

DCB Thunder, Fempac, Levant 1, Pacifier, IN.PACT 78% at 12m

Leipzig, Zeller Registry, IN.PACT 73-95% 12m

DES Zilver PTX, Real Zilver Ptx., Strides 85% at 12m

STENT grafts Viabahn 25cm 67% 12m

VIASTAR CS vs BMS 78 vs 53% 12m

Atherectomy alone and with DCB



LINC

Clinical Trials

Critical to the Research process
Well-controlled, Rigorously conducted

Limited by sample size
Patient selection criteria
Closer level of oversight
Training of clinical personnel

May not tell the full story about function in “the real world.”



LINC

Prospective Cohort

One Surgeon
Multiple hospitals
Enrolment 2016, Continuing
Research Officer Assessment
Data Tool: Pre, Intra and Post Procedure
FUP Discharge, 1m, 3m ,6m, 12m, and 6mthly
Antiplatelet 3m dual then single for life

Clinical Assessment, Duplex scan, ABI
Angiogram: if Clinical indicated

| | |
|------------|-------------|
| Rutherford | 2-6 |
| Lesions | Fem POP |
| | TASC A to D |
| | Restenosis |

No Clinical Events Committee
No Imaging Core Labs



L I N C

Outcomes

| | |
|---------|---|
| Patient | Mortality Early and late, Major Amputation, Procedure related |
| Lesion | Symptomatic : Clinical reintervention(CD TLR) Non Symptomatic: Primary Patency |

Major adverse event all-cause mortality, CD TVR, major target limb amputation, thrombosis at the target lesion site.

Primary safety composite endpoint was freedom from device- and procedure-related mortality through 30 days, and freedom from major target limb amputation and CD TVR within followup period



L I N C

Patient Demographic/comorbidities

| | (N=448) |
|--|--------------|
| Age (Y) | 63.1 ± 12.3 |
| Male Gender (%) | 62% |
| Diabetes (%) | 40% |
| Hypertension (%) | 75% |
| Hyperlipidemia (%) | 70% |
| Prior/ Current Smoker (%) | 65% |
| Coronary Artery Hx/Interv | 35% |
| Carotid Artery Hx/Interv.(%) | 25% |
| Renal Insufficiency (%) (cr>110umol/l) | 30% |
| Previous Peripheral Revasc. (%) | 55% |
| Tissue loss (%) | 15% |
| ABI | 0.609± 0.195 |



L I N C

Lesion Characteristic

| | |
|--------------------------------------|-----------------|
| Single | 84% |
| Multiple | 16% |
| <u>Lesion Type:</u> | |
| de novo | 74% |
| restenotic (no ISR) | 16% |
| ISR | 10% |
| Lesion Length | 13.97 ±12.57 cm |
| Total Occlusions | 19.8% |
| Severe Calcification | 20.1% |
| Diameter Stenosis (pre-treatment) | 80.8% ± 15.5 |
| No of balloon (1-8) | 2 ± 1.04 |

Procedure Outcomes

| | |
|-------------------|-------|
| Pre-dilatation | 100% |
| Device Success | 98% |
| Procedure Success | 100% |
| Clinical Success | 98.5% |
| Additional Stent | 25% |
| Distal embolus | 5.2% |
| Dissection | 25% |
| Rupture | 7% |

Followup Outcomes

Median 21m (1-37)

Fup 45pts (10.1%) lost

| Characteristic | Patient |
|------------------------------|---------|
| Clinically-driven TLR | 13.8% |
| Death (all-cause)(cardiac) | 6.9% |
| Major Target Limb Amputation | 1.1% |
| Thrombosis | 2.9% |
| Primary Safety Endpoint | 84% |
| Major Adverse Events | 24.7% |



LINC

Conclusions

Durable Safe and Effective IN.PACT Admiral in SFA/pop disease
in the Australian population in real world clinical practice

Problem single Arm Trial , no Independent monitor

Expensive treatment

Currently continuing to use DCB

Long term use depends on conclusion of current mortality controversy and PTX



LINC

**IN.PACT DCB results from
Real world experience of a
Single vascular surgeon
Prospective registry**

Kishore Sieunarine
Vascular Surgeon
Perth WA